


**AccelStim™**  
Bone Healing Therapy

**PhysioStim™**  
Bone Healing Therapy

## Fracture Management Solutions



**ORTHOFIX®**





## Accelerated Healing for Indicated Fresh Fractures

The AccelStim™ device provides a safe and effective nonsurgical treatment to improve nonunion fracture healing and accelerate the healing of indicated fresh fractures.\* The device stimulates the bone's natural healing process by low-intensity pulses of ultrasound (LIPUS) waves to the fracture site.<sup>1-4</sup>

### Proven Effective Therapy

- Accelerates Fracture Healing Recovery by 38%<sup>1</sup>
- Overall clinical success rate of 86% for nonunion fractures<sup>2</sup>
- 20 minutes daily treatment time

### Faster healing in both cortical and cancellous bone




### LIPUS has been proven to stimulate bone healing at the molecular, cellular, and tissue level.

#### Molecular:

LIPUS stimulates signaling pathways, leading to increased cell differentiation<sup>5,6</sup>

#### Cellular:

LIPUS stimulates bone cells to proliferate, differentiate, and mineralize<sup>7</sup>



#### Tissue:

LIPUS increases new bone formation<sup>8</sup>



Tibia/Fibula 89.6%<sup>9</sup>



Hand/Wrist 91.8%<sup>10</sup>




Scaphoid 92.2%<sup>11</sup>

**MOST RECENT**  
BONE GROWTH STIMULATION DEVICE  
ON THE MARKET

## For Patients Recovering from Nonunion Fractures

PhysioStim™ devices provide a safe and effective non-surgical treatment to improve nonunion fracture healing. These devices use a pulsed electromagnetic field (PEMF) signal to induce a low-level electrical field at the fracture site which stimulates bone healing.<sup>12,13</sup>


### Proven Effective Therapy



Healing Success in Common Nonunion Sites <sup>14</sup>	
Femur	84.2%
Fibula	91.4%
Metatarsal	90.9%
Tibia	89.0%
Ulna	96.1%
Radius	93.8%

### Anatomically Designed Models

- 360 degrees of PEMF treatment around the fracture site that evenly penetrates across tissue, bone and fixation<sup>15,16</sup>
- Single-piece, cordless design that allows for ease of placement and patient mobility
- Effective delivery of treatment when worn over clothing, casts and boots



Model 5313



Model 5314



Model 5302/5303



Model 5315



Model 5302/5303

### Commitment To Outcomes



STIM onTrack™  
Mobile App

PhysioStim devices are accompanied by the STIM onTrack mobile app. The app includes a first-to-market feature that:

- Enables physicians to remotely view patient adherence to their prescription
- Engages patients in their recovery process through treatment calendars, therapy reminders and educational resources

#1

ORTHOFIX PEMF BONE GROWTH THERAPY DEVICES

**PRESCRIBED**

BONE GROWTH STIMULATORS<sup>17,18</sup>



Orthofix Bone Growth Therapy devices are prescribed with a Guarantee Program which states that when prescribed for an approved FDA indication and when other eligibility requirements are met, radiographic progress will be shown in fracture healing or fusion healing, or the fee paid for the unit will be refunded to the payer(s) of record\*\* or, at the direction of the originally prescribing physician, a one-time replacement unit can be provided. This permits physicians to prescribe and insurance providers to approve our bone growth therapy devices with confidence, and most importantly, to assure our patients will have the maximum opportunity to heal.

**\*\*Subject to eligibility requirements.**

**\*Brief Prescribing Information:**

**The AccelStim™ device:** is indicated for the non-invasive treatment of established nonunions excluding skull and vertebra, and for accelerating the time to a healed fracture for fresh, closed, posteriorly displaced distal radius fractures and fresh, closed, or Grade I open tibial diaphysis fractures in skeletally mature adult individuals when these fractures are orthopedically managed by closed reduction and cast immobilization.

**The PhysioStim™ device:** is indicated for the treatment of an established nonunion acquired secondary to trauma, excluding vertebrae and all flat bones, where the width of the nonunion defect is less than one-half the width of the bone to be treated. A nonunion is considered to be established when the fracture site shows no visibly progressive signs of healing.

Use of this device is contraindicated where the individual has synovial pseudarthrosis. Demand type pacemaker operation may be adversely affected by exposure to pulsed electromagnetic fields. The safety and effectiveness of this device has not been established for individuals lacking skeletal maturity or individuals with a nonunion secondary to, or in connection with, a pathological condition. The safety of this device for use on patients who are pregnant or nursing has not been established. Rare instances of reversible minor discomfort have been reported.

Full prescribing information can be found in product labeling on our patient education website [BoneGrowthTherapy.com](http://BoneGrowthTherapy.com) or by calling Patient Services at 1-800-535-4492.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

**References:**

1. Kristiansen TK, Ryaby JP, McCabe J, Frey JJ, Roe LR. Accelerated healing of distal radial fractures with the use of specific, low-intensity ultrasound. *J Bone Joint Surg.* 1997;79-A(7):961-973. 2. Nolte PA, van der Krans A, Patka P, Janssen IMC, Ryaby JP, Albers GHR. Low-intensity pulsed ultrasound in the treatment of nonunions. *J Trauma.* 2001;51(4):693-703. 3. Heckman JD, Ryaby JP, McCabe J, Frey JJ, Kilcoyne RF. Acceleration of tibial fracture-healing by non-invasive, low-intensity pulsed ultrasound. *J Bone Joint Surg.* 1994;76-A(1):26-34. 4. PMA P210035. May 2022. 5. Parvizi J, Wu CC, Lewallen DG, Greenleaf JF, Bolander ME. Low-intensity ultrasound stimulates proteoglycan synthesis in rat chondrocytes by increasing aggrecan gene expression. *J Orthop Res: official publication of the Orthopaedic Research Society.* 1999;17(4):488e494. <https://doi.org/10.1002/jor.1100170405>. PubMed PMID: 10459753. 6. Azuma Y, Ito M, Harada Y, Takagi H, Ohta T, Jingushi S. Low-intensity pulsed ultrasound accelerates rat femoral fracture healing by acting on the various cellular reactions in the fracture callus. *J Bone Miner Res.* 2001;16(4):671-80. 7. Gurkan UA, Akkus O. The mechanical environment of bone marrow: a review. *Ann Biomed Eng.* 2008; 36(12): 1978e1991. <https://doi.org/10.1007/s10439-008-9577-x>. PubMed PMID: 18855142. 8. Wang, Y.; Peng, W.; Liu, X.; Zhu, M.; Sun, T.; Peng, Q.; Zeng, Y.; Peng, B.; Zhi, W.; Weng, J.; et al. Study of bilineage differentiation of human-bone-marrow-derived mesenchymal stem cells in oxidized sodium alginate/N-succinyl chitosan hydrogels and synergistic effects of RGD modification and low-intensity pulsed ultrasound. *Acta Biomater.* 2014, 10, 2518–2528. 9. Bioventus LLC. Tibia nonunion claims based on EXOGEN Registry. Data on file, RPT-000391. 10. Bioventus LLC. Hand/Wrist nonunion claims based on EXOGEN Registry. Data on file, RPT-0001411. 11. Bioventus LLC. Scaphoid nonunion claims based on EXOGEN Registry. Data on file, RPT-000398. 12. PMA P850007. February 1986. 13. Garland DE, Moses B, Salver W. Fracture healing: Long-term follow-up of fracture nonunions treated with PEMFs. *Contemp Orthop.* 1991;22(3):295-302. 14. Orthofix patient registry. PMA P850007/S20. Data on file. 15. Data on file. Field mapping analysis conducted by M. Zborowski, Ph.D., Cleveland Clinic. 16. Navarro, M., Michiardi, A., Castano, O., & Planell, J. (2008). Biomaterials in orthopaedics. *Journal of the Royal Society Interface*, 5(27), 1137-1158. 17. iData Research Inc., U.S. Market for Spinal Implants and VCF (iDATA\_USSP2021\_RPT), iData Research Inc (www.idataresearch.net) 2021. 18. iData Research Inc., U.S. Market for Orthopedic Trauma Devices (iDATA\_USTRA2021\_RMS), iData Research Inc (www.idataresearch.net) 2021.



STIM onTrack™  
Mobile App



Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

STIM onTrack is only available with PhysioStim™



[info@gobtwnow.com](mailto:info@gobtwnow.com)

479-271-2288

B5-2239 © Orthofix US LLC 01/2023

